

Bath Oil Taste

Bath for Oil Taste



Thermostated digital bath for OIL testing,
 controlled by microprocessor.

- Bath with electronic digital regulation, controlled by microprocessor.
- Aluminium block with 6 or 9 cavities, specially designed for the edible oil taste (standard no: L.248/67).
- Inner tank of stainless steel AISI 304 and metallic housing protected with Epoxi paint of great resistance to aggressive surroundings
- Digital Display LCD backlighted showing the selected and actual temperature, and timer.
- Membrane push-button control.
- Heating resistance of special alloy of stainless steel heat and corrosion resistance.
- Safety thermostat cut out indicator lamp. The unit will cut out if there is insufficient liquid or if the temperature exceeds the maximum temperature. DIN 12877 class 2.
- Indication on the display of high temperature on the tank, that provides protection against burn.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display.
- Memory of the last settings.

Technical specifications:

Model	BCA6E	BCA9E
Reference	10000-01133	10000-01135
Thermostatization		
Nº de cavities	6	9
T ^a range (°C)	T ^a amb+5 – 110	
Resolution (°C)	0,1	
Timer (min)	0-120	
Power (W)	600	
General Data		
Dimensions (mm)	270x380x160	
Weight (kg)	6,2	10,0
Protection IP	IP52	
Packaging cm	47x37x32	