

# Horizontal Maxi

## Digital Horizontal Shaker



Digital horizontal shaker, controlled by microprocessor.

- Horizontal (reciprocating) shaker, suitable for large size containers.
- Speed adjustable, controlled by microprocessor.
- Constant speed independently of the load.
- Smooth and progressive start and stop, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected and actual speed and timer (with audible alarm).
- Membrane push-button control.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display.
- Non-slipping surface, optional.
- Quiet operation. Suitable for continuous running.
- Packaging: 69x53x29cm

### Technical Specifications:

Model	HL30E
Reference	10000-01012
<b>Shaking</b>	
Max Load (Kg)	10
Power (W)	50
Type of Movement	Horizontal
Amplitude (mm)	30
Speed Range (rpm)	30-220
Resolution (rpm)	5
Timer (min)	0-120
<b>General Data</b>	
Plate Material	Non-slipping surface, optional
Useful surface (mm)	350x400
Dimensions (mm)	410x570x140
Weight (kg)	12,0
Protection IP	IP54

Accessories	Reference
<b>Universal Platform of 350 x 400 mm</b> To fix flask, precipitation vessels, bottles, grids of test tubes, etc	20000-00041
<b>Double Universal Platform of 350 x 400 mm</b> Same as the previous one, but with 2 shelves and 10 rollers	20000-00131
<b>Non-slipping Platform of 350 x 400 mm</b> Suitable for Petri plates, vessels, etc	20000-00084
<b>Double Non-slipping Platform of 350x400 mm</b> Same as the previous one, but with 2 shelves	20000-00157
<b>Platform for fastening Erlenmeyer flasks clips</b>	20000-00040
Clip for Erlenmeyer 50 ml. (up to 30u.)	20000-00046
Clip for Erlenmeyer 100 ml. (up to 20u.)	20000-00047
Clip for Erlenmeyer 250 ml. (up to 10u.)	20000-00048
Clip for Erlenmeyer 500 ml. (up to 8u.)	20000-00049
Clip for Erlenmeyer 1000 ml. (up to 5u.)	20000-00099
Clip for Erlenmeyer 2000 ml. (up to 4u)	20000-00209
<b>Magnetic platform for fastening Erlenmeyer flasks clips</b>	20000-00191
Magnetic clip for Erlenmeyer 50 ml. (up to 30u.)	20000-00192
Magnetic clip for Erlenmeyer 100 ml. (up to 20u.)	20000-00193
Magnetic clip for Erlenmeyer 250 ml. (up to 12u.)	20000-00194
Magnetic clip for Erlenmeyer 500 ml. (up to 9u.)	20000-00195
Magnetic clip for Erlenmeyer 1000 ml. (up to 6u.)	20000-00196
<b>Platform for Separating Funnels</b> For 5 separating funnels of 250 cc, 5 separating funnels of 500 cc. or 3 separating funnels of 1.000 cc	20000-00044

