

# MiniMagMix

## Magnetic Stirrer with Heating



Magnetic stirrer with heating, controlled by microprocessor.

- Magnetic stirrer with heating of reduced dimensions.
- LED's display with indication of the selected speed and temperature.
- Membrane push-button control.
- Appropriate for all kind of laboratories.
- Selection among 4 levels each one.
- Plate made of polished aluminium.
- Speed and temperature controlled by microprocessor.
- Packaging: 36x26x18

### Technical Specifications:

Model	MNG02E
Reference	10000-01033
<i>Magnetic Stirring</i>	
Max. Volume (L)	2
Power (W)	2,8
Speed Range (rpm)	300-1250
Resolution (rpm)	4 levels: 300-600-900-1250
<i>Heating</i>	
Heating Power (W)	250
Temperature Range (°C)	100-400
Resolution (°C)	4 levels: 100-200-300-400
<i>General Data</i>	
Plate Material	Polished aluminium
Plate Diameter (mm)	100x100
Dimensions (mm)	115x215x90
Weight (kg)	1,3
Protection IP	IP32

Accesorios	Referencia
<b>Set of 18 assorted magnets</b>	20000-00026
<b>Stirring magnet 8x30mm* (pack 5)</b>	20000-00008
* for other sizes, ask	

