

# MultiMix D

## Digital Magnetic Stirrer Multiplace without Heating



Digital magnetic stirrer multiplace without heating, controlled by microprocessor.



- Module consisting in 5 or 9 stirring places, in a single platform, with a single electrical connection. Each magnetic stirrer works independently.
- Adjustable stirring speed, controlled by microprocessor.
- Smooth and progressive start to maintain the magnetic coupling, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected speed.
- User friendly interface to adjust speed on each position.
- Membrane push-button control.
- The selected parameters are memorized, even if the instrument is switch off.
- Antiscratch surface.
- Packaging: 56x48x21cm

Technical Specifications:

Model	MM50E	MM90E
Reference	10000-01038	10000-01039
Nº places	5	9
<b>Magnetic Stirring</b>		
Max. Volume (L) (x plaza)	2	
Power (W)	9	16
Speed range (rpm)	200-1200	
Resolution (rpm)	100	
<b>General Data</b>		
Plate material	Stainless Steel AISI	
Plate diameter (mm) (x plaza)	110	
Dimensions (mm)	365x445x60	
Weight (kg)	3,5	3,8
Protection IP	IP53	

Option of Stirrer + Tank + Immersion Heater	Reference
<i>Bath OvanTherm MultiMix BHM9E (9 places)</i>	10000-01128
<i>Bath OvanTherm MultiMix BHM5E (5 places)</i>	10000-01130



Accessories	Reference
<b>Set of 18 assorted magnets</b>	20000-00026
<b>Stirring magnet 8x30mm * (pack 5)</b> *for other sizes, ask	20000-00008

