

MultiMix Heat

Magnetic Stirrer Multiplace with Heating



Magnetic stirrer multiplace with heating, controlled by microprocessor.

- Module consisting of 3 stirring units in a single platform, with a single electrical connection. Each magnetic stirrer works independently.
- Speed and temperature adjustable, controlled by microprocessor.
- Selection among 10 different levels for temperature as well as for speed.
- Display of LED's with indication of the selected values.
- Membrane push-button control.
- Removable plates to be able to replace itself in case of deterioration.
- Packaging: 63x32x22cm

Technical Specifications:

Model	MMH3E
Reference	10000-01022
Nº places	3
Magnetic Stirring	
Max. Volume (L) (x place)	5
Power (W) (x place)	5
Speed Range (rpm)	300-1500
Resolution (rpm)	10 levels
Heating	
Heating Power (W) (x place)	500
T ^a range (°C)	37-400
Resolution (°C)	10 levels
Connection to ext. T ^a sensor	Contact Thermometer or Electronic Thermometer (need to be ordered)
General Data	
Plate material	Aluminium alloy
Plate diameter (mm)	154 (x place)
Dimensions (mm)	525x205x120
Weight (kg)	7
Protection IP	IP32

Accessories	Reference
Contact Thermometer 100°C 15cm rod	20000-00185
Contact Thermometer 200°C 15cm rod	20000-00186
Contact Thermometer 300°C 15cm rod	20000-00187
Cable for thermometer DIN connection	20000-00088
Support Rod Includes 2 fastening wing nuts to fix it. 450 mm	20000-00036
Fastening clamp for Pt100	20000-00073
Set of 18 assorted magnets	20000-00026
Stirring magnet 8x30mm* (pack 5) *for other sizes, ask	20000-00008
Protection ring stainless steel	20000-00028

