

RollerMix D

Digital Rotary and Oscillating Mixer with Rollers



Digital rotary and oscillating mixer with rollers, controlled by microprocessor.

- Mixer with simultaneous double rotative-oscillating movement.
- For test tubes, bottles and closed containers.
- Speed adjustable, controlled by microprocessor .
- Five agitation rollers forming four available rails.
- Digital display LCD backlighted indicating the selected and actual speed
- and timer (with audible alarm).
- Membrane push-button control
- Geared synchronization with no wear.
- Brushless motor with unblocking system by clutch.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display
- Quiet operation, suitable for continuous running.
- Packaging: 70x31x22cm

Technical Specifications:

Model	RM120-DE
Reference	10000-01071
Mixing	
Max.Load(kg)	1,5
Power (W)	4
Type of Movement	Rotation – Oscillation*
Nº rollers	5
Inclination angle *	+/-4º
Speed range (rpm)	30-150
Resolution (rpm)	5
Timer (min)	0-120
General Data	
Length of rollers (mm)	320
Dimensions (mm)	550x190x80
Weight (kg)	3,5
Protection IP	IP54

*can be ordered without oscillation movement