

# Swing

## Digital Rocking Mixer



Digital rocking shaker, controlled by microprocessor.

- Speed adjustable, controlled by microprocessor.
- Constant speed independently of the load.
- Smooth and progressive start and stop, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected and actual speed and timer (with audible alarm).
- Membrane push-button control.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display.
- Non-slipping surface.
- Quiet operation. Suitable for continuous running.
- Packaging: 52x46x25cm

Technical Specifications:

| Model                | SW8E                 |
|----------------------|----------------------|
| Reference            | 10000-01015          |
| <i>Shaking</i>       |                      |
| Max. Load (Kg)       | 5                    |
| Power (W)            | 50                   |
| Type of movement     | Rocking              |
| Angle of slope       | +/-8°                |
| Speed Range (osc.pm) | 5-50                 |
| Resolution (osc.pm)  | 1                    |
| Timer (min)          | 0-120                |
| <i>General Data</i>  |                      |
| Plate Material       | Non-slipping surface |
| Useful surface (mm)  | 290x260              |
| Dimensions (mm)      | 340x400x125          |
| Weight (kg)          | 4                    |
| Protection IP        | IP54                 |

| Accessories   | Reference   |
|---|-------------|
| <b>Mat with 12 channels, 260 x 290 mm</b><br>For tubes of different diameters | 20000-00056 |