

VibraMix

Vortex Mixer



Vortex mixer shaker.

- Universal mixer with orbital motion and vibration effect thanks to the small orbit.
- Appropriate for test tubes, vessels, etc.
- Adjustable and variable speed with external potentiometer, controlled by microprocessor.
- Constant speed independently of the load.
- Two modes of operation:
 - continuous
 - touch start running.
- Wide header that allows to shake tubes from 6 to 30 mm Ø.
- Variety of accessories for different applications.
- Automatic restart after unexpected shut down.
- For its hardness and reliability, it doesn't need suction pads, and offers a major reliability.
- Packaging: 20x26x18cm

Technical Specifications:

Model	VM3E
Reference	10000-01006
Shaking	
Power (W)	12
Type of movement	Orbital
Amplitude (mm)	3
Speed range (rpm)	0-2600
General Data	
Dimensions (mm)	125x170x110
Weight (kg)	4
Protection IP	IP54

Accessories	Reference
Platform Ø 90 mm To shake vessels or several tubes at the same time	20000-00042
Header Ø 60 mm For 10 Eppendorf tubes Ø 10 x 50 mm	20000-00083
Support shaft for tubes * length 60 mm	20000-00086
Support rod * Includes 2 fastening wing nuts to fix it. 450 mm. length.	20000-00036



* Together with the two previous references enables assemblies for continuous operation