

IV OVANGEL

VERTICAL LAMINAR FLOW CABINET



IV Ovangel is a laminar flow cabinet with a clean working area of Class ISO 5 (ISO 14644).Working function by sucking air through a G-4 prefilter and blowing it through a Hepa H-14 filter .

- Vertical Laminar flow: ISO 5
- Hepa filter with Gel sealant fitting
- Low noise thanks to DC motor
- Modular concept
- CPT digital control

Technical Specifications:

Model	IV Ovangel 9	IV Ovangel 12
Reference	50000-00121	50000-00122
Overall Dimensions (mm)	930 x 1190 x 610	1230 x 1190 x 610
Useful Dimensions (mm)	900 x 660 x590	1190 x 660 x 590
Weight (Kg)	90	100
Nominal Power (W)	300	300
Electr.Connection (V – Hz)	230 -115 V-50/60Hz, single phase	
Serial Equipment		
Prefilter	W Type	
Filter	Minipleat: 99,995% MPPS – EN- 1822μ	
Light	Fluorescent: 400 – 600 lux.	
Fan	ECM	

Options	Reference
<i>Epoxi painted base for IV OvanGel</i>	50000-01118
<i>Other options available .- under request</i>	