

VERTIOVAN MV 9

VERTICAL FLOW CABINET



Air is taken from the top through a G4 pleated prefilter. The air is sucked by a very silent engine and then blown into a Minipleat Absolute Filter HEPA directly to the glass work area.

- Security glass work area
- Plexi glass top flip window
- Class 100
- Blue Gel conception
- Analogue control

Technical Specifications:

Model	Vertiovan MV 9
Reference	50000-00120
Overall Dimensions (mm)	925 x 1030 x 600
Useful Dimensions (mm)	905 x 700 x 592
Weight (Kg)	30
Electr.Connection (V – Hz)	220V-50Hz, single phase
Serial Equipment	
Prefilter	W Type
Filter	Minipleat: 99,999% 0,3 μ.
Light	Fluorescent: 400 – 600 lux.
Fan	Speed variation

Options	Reference
<i>UV Lamp Vertiovan MV 9</i>	<i>50000-01108</i>
<i>Resin counter top with steel base</i>	<i>50000-01109</i>
<i>Other options available .- under request</i>	