

Ultrasonic Bath

Thermostated ultrasonic cleaner 2,5 l



Thermostatic Ultrasonic Bath

- Thermostated Ultrasound bath.
- S. Steel inner tank
- Front part anti spill
- Temperature controlled by microprocessor.
- Integrated control display
- Drain valve (except 0,7 and 2,5 cleaner).

Technical data:

Reference	40000-00300	40000-00026*	40000-00301	40000-00302	40000-00303	40000-00304	40000-00305
Vol. Máx. (L)**	1,5	2,5	3	6	10	14	28
Tª (°C)	Tªamb -80	Tamb -65°C	Tªamb - 80				
Frequency kHz	40	42	40				
Precision (°C)	0,1	0,1	0,1				
Ultrasound power(W)	50	170	100	200	300	500	600
Heating power (W)	100	100	100	150	200	300	750
General data							
Inner tank	St.Steel AISI 304	St.Steel AISI 304	St.Steel AISI 304				
Useful dimensions (mm)	140 x 130 x 100	264x164x80	230 x 130 x 100	300 x 150 x 150	295 x 235 x 150	325 x 295 x 150	495 x 295 x 200
Weight (kg)	2,3	2,5	3,7	5,7	9	14	18

Referencia	40000-00306	40000-00307
Vol. Máx. (L)**	35	50
Tª (°C)	Tªamb -80	
Frequency kHz	40	
Precision (°C)	0,1	
Ultrasound power(W)	700	800
Heating power (W)	1000	1000
General data		
Inner tank	St.Steel AISI 304	
Useful dimensions (mm)	395 x 295 x 300	600 x 300 x 300
Weight (kg)	25	40

* 2,5 l plastic external housing

** Different volumes under request

- Timer: 0- 30 min (with a beep sound alert)
- Lid included
- Basket included
- Thermostat from Ambient temp to 80°C