

# OPTIOVANGEL

## HORIZONTAL LAMINAR FLOW CABINET



Optiovangel provides a clean air working area of Class ISO 5 (ISO 14644). Air is constantly blown horizontally through an absolute HEPA filter into the working area, protecting the product inside the cabinet from existing particles in the air, preventing cross-contamination.

- Horizontal Laminar flow: ISO 5
- HEPA filter H-14
- Low noise thanks to DC motor
- CPT digital control

**Technical Specifications:**

<b>Model</b>	<b>OptioOvangel 12</b>	<b>OptioOvangel 18</b>
<b>Reference</b>	50000-00125	50000-00126
<b>Overall Dimensions (mm)</b>	1243 x 800 x 850	1875 x 800 x 850
<b>Useful Dimensions (mm)</b>	1212 x 620 x 540	1815 x 620 x 540
<b>Weight (Kg)</b>	100	150
<b>Nominal Power (W)</b>	150	150
<b>Electr. Connection (V – Hz)</b>	220V-50Hz, single phase	
<b>Serial Equipment</b>		
<b>Prefilter</b>	2 prefilters W Type G-3	
<b>Filter</b>	Minipleat HEPA H-14	
<b>Light</b>	Fluorescent: 400 – 600 lux.	
<b>Fan</b>	ECM	

<b>Options</b>	<b>Reference</b>
<i>Epoxi painted base for OptiOvanGel</i>	<b>50000-01110</b>
<i>Other options available .- under request</i>	