

Jar Test

Digital Flocculator



Digital flocculator controlled by microprocessor.



- Digital flocculator for water analysis of 6 or 4 places.
- Suitable for long vessels of 1L.
- Speed controlled by microprocessor.
- Smooth and progressive start and stop.
- Digital display LCD backlighted indicating the selected and actual speed and timer (audible alarm).
- Membrane push-button control.
- Memory for 2 different speed values, changing between each one just with one click..
- Stainless steel stirring rods of Ø 6 mm. with adjustable height.
- Backlight screen that could be placed in 2 positions: lower or rear, without dismantling.
- Handles for easy transport.
- Quiet operation, suitable for continuous running.
- Packaging: Tribox: 84x50x32cm

Technical specifications:

Model	JT40E	JT60E
Reference	10000-01016	10000-01017
Mixing		
No places	4	6
Power (W)	50	
Speed range (rpm)	20-250	
Resolution (rpm)	5	
Accuracy (rpm)	5	
Timer (min)	0-240	
General Data		
Dimensions (mm)	800x200x410	
Weight (kg)	14,0	15,0
Protection IP	IP34	

Accessories	Reference
Injector of coagulants for reagent (4 and 6 places)	20000-00079