

Oven van N

Natural convection Ovens

Digital Oven, controlled by microprocessor



- Electronic controller at the instrument lower part to avoid overheating
- Natural convection oven to avoid inner sample drying
- Steel Inner chamber
- Removable shelves
- Epoxi paint external housing

Specifications:

Model	Ovenvan N 24	Ovenvan N 67
Referencia	59000-00600	59000-00601
Thermostatization		
Vol. (l)	24	67
Temp range(°C)	T amb +10 - 200	T amb +10 - 350
T estability(± °C)	±1	±1
T Homogeneity at 300°C	± 5	± 14
Voltage (V)	230	230
Potwer (W)	1000	2000
General Data		
Inner chamber material	Steel	
Shelves number	2	3
Max shelves	2	7
Max weight by shelf (Kg)	10	10
Useful dimensions. (mm) (Wide x depth x height)	300x380x200	390x445x390
External dimensions (mm) (Wide x depth x height)	400x515x410	675x625x600
Net Weight (kg)	24	37

Accessories	
59000-01500	Buzzer
59000-01501	Timer
59000-01502	OTP- Temp overprotection
59000-01504	Programmable controller
59000-01505	Interface RS 232 or 485
59000-01506	Pc connection software
59000-01507	Sample metallic box
59000-01508	Inlet gas box .