

Jar Test P

Digital Portable Flocculator



Digital portable flocculator controlled by microprocessor.

- Digital portable flocculator for water analysis of 4 places.
- Suitable for long vessels of 1L.
- 2 working modes:
 - Free and adjustable stirring speed
 - **Pre-programmed under ASTM D-2035 conditions**
- Smooth and progressive start and stop.
- Digital display LCD backlighted indicating the selected and actual speed and timer (audible alarm).
- Membrane push-button control.
- Memory for 2 different speed values, changing between each one just with one click..
- Stainless steel stirring rods of Ø 6 mm. with adjustable height.
- Light screen at each position.
- Handle for easy transport.
- Double power supply: 230V AC and 12V DC, that allows to use it with a car battery.
- Includes a cable for the car lighter connector.
- Packaging: 66x45x43cm

Technical specifications:

Model	JT4E
Reference	10000-00097
Mixing	
No places	4
Power (W)	50
Speed range (rpm)	20-250
Resolution (rpm)	5
Accuracy (rpm)	±1
Timer (min)	0-240
General Data	
Dimensions (mm)	250x250x290
Weight (kg)	6,0
Protection IP	IP34