

MicroMagMix

Digital Magnetic Stirrer with Heating



Digital magnetic stirrer with heating, controlled by microprocessor.



- Appropriate for all kind of laboratories, thanks to its robustness and its ergonomic and modern design.
- Speed and temperature controlled by microprocessor, by means of an integrated Pt100 sensor.
- ARA V2 fuzzy logic algorithm for high precision temperature control, to be used with external sensor, without surpassing reference temperature (even at 30° C).
- Constant speed independently that the conditions of the sample vary.
- Smooth and progressive start to maintain the magnetic coupling, independently of the selected speed.
- Digital display LCD backlighted indicating the selected and actual speed and temperature.
- Indicator on the display of high temperature on the plate, that provides protection against burns.
- Membrane push-button control.
- Automatic restart after unexpected shut down, keeping the same parameters as before the incident. An incident sign appears on the display.
- Removable plate, to replace it in case of deterioration.
- Packaging: 36x26x18cm

Technical Specifications:

Model	MCG05E	MCG15E
Reference	10000-00092	10000-00093
<i>Magnetic Stirring</i>		
Max. Volume (L)	5	15
Power (W)	5	12
Speed Range (rpm)	100-1250	200-1500
Resolution (rpm)	25	
<i>Heating</i>		
Heating power (W)	500	
Temperature Range (°C)	20-300	
Resolution (°C)	0,5	
Accuracy (°C)	+/- 1	
Max. Viscosity (cSt)	10	40
Connection to extern T ³ sensor	Pt100	
<i>General Data</i>		
Plate Material	Polished Aluminium	Stainless Steel AISI 304
Plate Diameter (mm)	154	
Dimensions (mm)	210x260x130	
Weight	4	4,2
Protection IP	IP52	

Available Kits		Reference
<i>Kit MCG05E</i>	Includes MicroMagMix, extern Pt100 sensor, support rod and fastening clamp	10000-00098
<i>Kit MCG15E</i>		10000-00099

Accessories	Referencia
External Pt100 Sensor Allows temperature control inside the preparation	20000-00034
Support Rod Includes 2 fastening wing nuts to fix it. 450 mm.	20000-00036
Fastening clamp for Pt100	20000-00073
Caja surtida 18 imanes	20000-00026
Set of 18 assorted magnets* (pack 5) *for other sizes, ask	20000-00008
Kjeldhal adapter 250cc (only for MCG15E)	20000-00032
Kjeldhal adapter 500cc (only for MCG15E)	20000-00033
Protection ring stainless steel	20000-00028

