

MicroMix

Digital Magnetic Stirrer without heating



Digital magnetic stirrer without heating, controlled by microprocessor.

- Appropriate for all kind of laboratories, thanks to its robustness and its ergonomic and modern design.
- Adjustable stirring speed, controlled by microprocessor.
- Constant speed independently that the conditions of the sample vary.
- Smooth and progressive start to maintain the magnetic coupling, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected and actual speed and timer (audible alarm).
- Membrane push-button control.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display.
- Packaging: 36x26x18cm

Technical Specifications:

| Model | MC08E | MC20E |
|--------------------------|--------------------------|-------------|
| Reference | 10000-00094 | 10000-00095 |
| Magnetic Stirring | | |
| Max. Volume (L) | 8 | 20 |
| Power (W) | 12 | 30 |
| Speed range (rpm) | 200-1500 | |
| Resolution (rpm) | 20 | |
| Timer (min) | 0-120 | |
| Max. Viscosity (cSt) | 10 | 40 |
| General Data | | |
| Plate material | Stainless Steel AISI 304 | |
| Plate diameter (mm) | 150 | |
| Dimensions (mm) | 210x260x110 | |
| Weight (kg) | 3 | 3,2 |
| Protection IP | IP52 | |

| Accessories | Reference |
|--|-------------|
| Support Rod Includes 2 fastening wing nuts to fix it. 450 mm | 20000-00036 |
| Set of 18 assorted magnets | 20000-00026 |
| Stirring magnet 8x30mm (pack 5) * for other sizes, ask | 20000-00008 |

