

MicroMix Ultra

Digital Magnetic Stirrer without Heating



Digital magnetic stirrer without heating, controlled by microprocessor.

- Adjustable stirring speed, controlled by microprocessor.
- Suitable for cinematic viscosity up to 100cSt.
- Constant speed independently that the conditions of the sample vary.
- Smooth and progressive start to maintain the magnetic coupling, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected and actual speed and timer (audible alarm).
- Membrane push-button control.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display.
- Very quiet, for continuous operation.
- Antiscratch surface of stirring.
- Packaging: 69x53x29cm

Technical specifications:

Model	MC50UE
Reference	10000-01052
Magnetic Stirring	
Max. Volume (L)	50
Power (W)	50
Speed range (rpm)	50-500
Resolution (rpm)	5
Timer (min)	0-120
General Data	
Plate material	Stainless steel AISI 304
Plate diameter (mm)	350x400
Dimensions (mm)	410x570x140
Weight	10
Protection IP	IP53

Accessories	Reference
Set of 18 assorted magnets	20000-00026
Stirring magnet 16x102 mm* * for other sizes, ask	20000-00000

