

MiniMix

Magnetic Stirrer without heating



Magnetic Stirrer without heating, controlled by microprocessor.

- Stirrer of reduced dimensions.
- Adjustable stirring speed, controlled by microprocessor.
- Membrane push-button control.
- Appropriated for all kind of laboratories.
- LED's display with indication of the selected speed in 4 levels.
- Antiscratch surface of stirring.
- Packaging: 36x26x18cm

Technical specifications:

Model	MN02E
Reference	10000-01036
Magnetic Stirring	
Max. Volume (L)	2
Power (W)	2,8
Speed range (rpm)	300-1250
Resolution (rpm)	4 levels: 300-600-900-1250
General Data	
Plate material	Aluminium alloy
Plate diameter (mm)	110
Dimensions (mm)	115x200x50
Weight	0,6
Protection IP	IP43

Accessories	Reference
Set of 18 assorted magnets	20000-00026
Stirring magnet 8x30mm * (pack 5) *for other sizes, ask	20000-00008

