

MultiMix Heat D

Digital Magnetic Stirrer Multiplace with Heating



Digital magnetic stirrer multiplace with heating, controlled by microprocessor.



- Module consisting in 5 or 9 stirring/heating places, in a single platform, with a single electrical connection. Each magnetic stirrer works independently.
- Adjustable speed and temperature, controlled by microprocessor.
- Smooth and progressive start to maintain the magnetic coupling, independently of the selected speed.
- Digital Display LCD backlighted indicating the selected speed and temperature of each place.
- User friendly interface to adjust speed on each position.
- Membrane push-button control.
- The selected parameters are memorized, even if the instrument is switch off.
- Hot plates made of polished aluminium.
- Packaging: 56x48x21cm
- 5p= 1250w
- 9p= 2300w

Technical Specifications:

| Model | MMH50E | MMH90E |
|-------------------------------|--------------------|-------------|
| Reference | 10000-01037 | 10000-01032 |
| Nº places | 5 | 9 |
| <i>Magnetic Stirring</i> | | |
| Max. Volume (L) (x place) | 2 | |
| Power (W) (x place) | 9 | 16 |
| Speed range (rpm) | 200-1200 | |
| Resolution (rpm) | 100 | |
| <i>Heating</i> | | |
| Heating power (W) (x place) | 250 | |
| Temperature range (°C) | 30-400 | |
| Resolution (°C) | 10 | |
| <i>General Data</i> | | |
| Plate material | Polished aluminium | |
| Plate diameter (mm) (x place) | 100x100 | |
| Dimensions (mm) | 360x445x90 | |
| Weight (Kg) | 6,0 | 7,0 |
| Protection IP | IP53 | |

| Accessories | Reference |
|--|-------------|
| Set of 18 assorted magnets | 20000-00026 |
| Stirring magnet 8x30mm* (pack 5) *for other sizes, ask | 20000-00008 |

