

Noria

Rotary mixer



Rotary mixer.

- For rotary motion with various tubes from 6 to 35 mm diameter.
- Adaptable to any tube dimension depending on the needs.
- Fix rotation speed
- Allows place and removal of the tubes without stops the machine.
- Completely adjustable slope from 0° to 90°.
- Security clutch allows to stop rotary motion with the hand without problems.
- Antiscratching surface with chemical resistant coating
- Suitable for continuous running.
- Packaging: 47x37x32cm

Technical Specifications:

Model	N12E
Reference	10000-01013
Mixing	
Type of movement	Rotary
Angle of slope	0°-90°
Speed Range (rpm)	12 (fix)
General Data	
Dimensions (mm)	320x260x200
Weight (kg)	4,7
Protection IP	IP54

Accessories	Reference
Platform with clamps for:	
25 tubes Ø 6-7 mm	20000-00093
25 tubes Ø 8-9 mm	20000-00094
25 tubes Ø 10-11 mm Eppendorf	20000-00095
20 tubes Ø 12-13 mm	20000-00133
20 tubes Ø 14-15 mm	20000-00134
20 tubes Ø 16-17 mm	20000-00135
19 tubes Ø 18-19 mm	20000-00136
16 tubes Ø 20-21 mm	20000-00137
14 tubes Ø 22-23 mm	20000-00138
14 tubes Ø 24-25 mm	20000-00139
14 tubes Ø 27-30 mm	20000-00140
10 tubes Ø 32-35 mm	20000-00141
Racks for:	
60 culture tubes Ø 12-13,5 mm	20000-00051
36 centrifuge tubes 15ml Ø 16-17,5mm	20000-00052
20 centrifuge tubes 50ml Ø 29-31mm	20000-00053
60 tubes Eppendorf de 0,5 a 2 ml	20000-00096

