

Noria RD

Digital Rotary mixer



Digital Rotary mixer, controlled by microprocessor.

- For rotary motion with various tubes from 6 to 35 mm diameter.
- Adaptable to any tube dimension depending on the needs.
- Rotation speed adjustable, controlled by microprocessor.
- LCD backlight display indicating the selected and actual speed and timer. Audible alarm.
- Membrane push-button control.
- Allows place and removal of the tubes without stops the machine.
- Completely adjustable slope from 0° to 90°.
- Smooth and progressive start and stop.
- Timer from 0 to 120 min.
- Security clutch allows to stop rotary motion with the hand without problems.
- Anti-scratching surface with chemical resistant coating.
- Quiet operation. Suitable for continuous running.

Technical Specifications:

Model	NR100DE
Reference	30000-00084
Mixing	
Power (W)	50
Type of movement	Rotary
Angle of slope	0°-90°
Speed Range (rpm)	10-100
General Data	
Dimensions (mm)	320x260x290
Weight (kg)	6,5
Protection IP	IP54

Accessories	Reference
Platform with clamps for:	
25 tubes Ø 6-7 mm	20000-00093
25 tubes Ø 8-9 mm	20000-00094
25 tubes Ø 10-11 mm Eppendorf	20000-00095
20 tubes Ø 12-13 mm	20000-00133
20 tubes Ø 14-15 mm	20000-00134
20 tubes Ø 16-17 mm	20000-00135
19 tubes Ø 18-19 mm	20000-00136
16 tubes Ø 20-21 mm	20000-00137
14 tubes Ø 22-23 mm	20000-00138
14 tubes Ø 24-25 mm	20000-00139
14 tubes Ø 27-30 mm	20000-00140
10 tubes Ø 32-35 mm	20000-00141
Racks for:	
60 culture tubes Ø 12-13,5 mm	20000-00051
36 centrifuge tubes 15ml Ø 16-17,5mm	20000-00052
20 centrifuge tubes 50ml Ø 29-31mm	20000-00053
60 tubes Eppendorf de 0,5 a 2 ml	20000-00096

