

RollerMix

Rotary and Oscillating Mixer with Rollers



Rotary and oscillating mixer with rollers.

- Mixer with simultaneous double rotative-oscillating movement.
- For test tubes, bottles and closed containers.
- Fix rotary speed.
- Five agitation rollers forming four available rails.
- Geared synchronisation with no wear.
- Brushless motor with unblocking system by clutch.
- Quiet operation, suitable for continuous running.
- Packaging: 58x32x22cm

Technical Specifications:

Model	RM20E
Reference	10000-01009
Mixing	
Max. Load (Kg)	2
Power (W)	4
Type of Movement	Rotation – Oscillation*
Nº rollers	5
Inclination angle *	+/-4º
Speed range (rpm)	20 (fix)
General Data	
Length of rollers (mm)	320
Dimensions (mm)	450x190x80
Weight (kg)	3,0
Protection IP	IP54

*can be ordered without oscillation movement