

Ultrasonic Bath

Ultrasonic Baths



Industrial series



Digital Ultrasonic Lab Bath

- Thermostatic Ultrasonic Baths.
- Inner Stainless Steel Tank
- Micropocessor control.
- Integrated LCD control display
- Exhaust valve and carry handles.

Ultrasonic Baths

Specifications:

Reference	40000-00300	40000-00301	40000-00302	40000-00303	40000-00304	40000-00305	40000-00306	40000-00307
Thermostatization								
Max. Volum (L) **	2	3	6	10	20	30	40	50
Temp range (°C)	Amb T - 80							
Freq. kHz	40							
Precision (°C)	0,1							
Ultrasound Powers(W)	50	100	200	300	500	600	600	800
Heating Power (W)	100	100	150	200	300	750	750	1000
Global Data								
Inner tank material	SS AISI 304							
Useful dimensions (mm)	150 x 135 x 100	240 x 130 x 100	300 x 150 x 150	300 x 240 x 150	495 x 295 x 150	495 x 295 x 200	500 x 350 x 250	600 x 300 x 300
Weight(kg)	2,3	3,7	5,7	9	14	16	40	43