

OVANHOOD MAV 9

CHEMICAL HOOD WITH CHARCOAL FILTER Ductless Fume Hood



Filtration is performed thanks to activated active carbon filter, responsible to absorb all gases and vapours present in the air flow going through.
Air is taken in the work volume thanks to the negative pressure sucked by a very silent motor/impeller, once is filtered the air is blown out on the top of the unit back to the room

- Security glass
- Plexi glass top slip Windows
- No duct
- Integrated prefilter
- Portable
- Silent

Technical Specifications:

Model	OvanHood 9
Reference	50000-00131
Overall Dimensions (mm)	925 x 1030 x 600
Useful Dimensions (mm)	905 x 700 x 592
Electr.Connection (V – Hz)	220V-50Hz, single phase
Serial Equipment	
Carbon Active Filter	NC35 type
Prefilter	85% Gravi Integrated
Light	White fluorescent tubes : 400 – 600 lux.
Fan	adjustable

Options	Reference
<i>Bactericide UV Lamp</i>	<i>50000-01115</i>
<i>Resin counter top with steel base</i>	<i>50000-01116</i>
<i>Charcoal filter: CALK 3/9</i>	<i>50000-01117</i>
<i>Other options available .- under request</i>	