

PPSOVAN

SAFETY WEIGHTING WORK STATION



PPSOvan is a laminar air flow Classe 100-ISO5, blows air through an absolute filter in the work space. This air flow is drawn through lateral grids going to a second Hepa filter beneath the work surface

- Laminar flow: ISO 5
- 2 x Hepa filter H-14
- Digital control
- Stainless Steel structure
- Automatic air flow regulation

Technical Specifications:

Model	PPSOvan 9	PPSOvan 12	PPSOvan 16	PPSOvan 18
Reference	50000-00127	50000-00128	50000-00129	50000-00130
Overall Dimensions (mm)	978x2260x950	1283x2260x950	1588x2260x950	1893x2260x950
Useful Dimensions (mm)	921 x650 x640	1226 x650 x640	1531 x650 x640	1836921 x650 x640

Options	Reference
<i>UV Lamp for PPOvan 9 and 12</i>	50000-01111
<i>UV Lamp for PPOvan 16 and 18</i>	50000-01112
<i>Gas inlet for PPOvan</i>	50000-01113
<i>Antivibration marble</i>	50000-01114
<i>Other options available .- under request</i>	