


Slide Drying

Analog and Digital Slide Dry



Slide Drying Units used in the preparation of microscope slides at preparation step.

Both units, analog and Digital have up to 52 slides capacity, and slides can be arranged in different orientations, hold or bending

Model	SD52-A	SD52-D
Reference	30000-00260	30000-00261
Heating		
Display	Analogic	Digital
Control	Analog Thermostat	Microprocessor
Nº Slides (hold/bend)	24/52	24/52
Heating power (W)	150	150
T ^a range (°C)	ambT ^a +5 - 90	
Resolution (°C)	-	0,5°C
General Data		
Useful surface (cm)	38 x 18	
Power Supply	230V 50Hz **	
Weight (kg)	2,2	2,3
Accessories Slide Drying		Reference
Rod		20000-0290