


Slide Drying

Analog and Digital Slide Dry



Slide Drying Units used in the preparation of microscope slides at preparation step.

Both units, analog and Digital have up to 50 slides capacity, and slides can be arranged in different orientations

Model	SD50-A	SD50-D
Reference	10000-XXXX	10000-XXXX
Heating		
Display	Analogic	Digital
Control	Analog Thermostat	Microprocessor
N° Slides	50	50
Heating power (W)	350	350
T ^e range (°C)	ambT ^a +5 - 110	
Resolution (°C)	-	0,1°C
General Data		
Power Supply	230V 50Hz **	
Weight (kg)	2,2	2,3