

Histo - Bath

Digital



TISSUE FLOTATION BATH FOR PARAFFIN SECTIONS

Black cuvette for better visualization of the cuts/sections

Tempered edge for sample drying and slide

Model	HBD
Reference	30000-10005
Heating	
Display	Digital
Control	Microprocessor
Power (W)	325
T ^a range (°C)	T amb+5 - 90
Resolution (°C)	0,5
General Data	
Measures (mm)	355 x 350 x 110
Fuse	3,15A
Power Supply	230V / 110V*
Weight (kg)	4,4