

# Jar Test

## Digital Flocculator



Digital flocculator controlled by microprocessor.



- Digital flocculator for water analysis of 6 or 4 places.
- Suitable for long vessels of 1L.
- 2 working modes:
  - Free and adjustable stirring speed
  - **Pre-programmed under ASTM D-2035 conditions**
- Smooth and progressive start and stop.
- Digital display LCD backlighted indicating the selected and actual speed and timer (audible alarm).
- Membrane push-button control.
- Memory for 2 different speed values, changing between each one just with one click..
- Stainless steel stirring rods of Ø 6 mm. with adjustable height.
- Backlight screen that could be placed in 2 positions: lower or rear, without dismantling.
- Handles for easy transport.
- Quiet operation, suitable for continuous running.
- Packaging: Tribox: 84x50x32cm

Technical specifications:

| Model               | JT40E       | JT60E       |
|---------------------|-------------|-------------|
| Reference           | 10000-01016 | 10000-01017 |
| <b>Mixing</b>       |             |             |
| No places           | 4           | 6           |
| Power (W)           | 50          |             |
| Speed range (rpm)   | 20-250      |             |
| Resolution (rpm)    | 5           |             |
| Accuracy (rpm)      | ±1          |             |
| Timer (min)         | 0-240       |             |
| <b>General Data</b> |             |             |
| Dimensions (mm)     | 800x200x410 |             |
| Weight (kg)         | 14,0        | 15,0        |
| Protection IP       | IP34        |             |

| Accessories  | Reference   |
|--|-------------|
| <b>Injector of coagulants for reagent</b><br>(4 and 6 places ) | 20000-00079 |