


Horizontal Platelet Shaker

Horizontal Digital Shaker



Digital Horizontal Shaker

- Horizontal Shaker
- Shaking speed fixed, controlled by microprocessor.
- Smooth start and stop .
- Friendly backlight LCD
- Membrane buttons at the keyboard.
- After power short cut unit begin working again **.
- Optional: stop/start delay button
- Included 9 shelves platform (560 x 400mm)

Technical Data

	Horizontal HBS
Model	HBS36RA
Reference	10000-01158

SHAKING	
Movement	HORIZONTAL
Amplitude (mm)	36
Max load (kg)	25
Max Output (W/N-m)	50/1
Speed range (rpm)	50-100
Resolution (rpm)	5
Timer Auto start delay (min)	optional
Timeer (min)	0-120*

General Data	
Volt. (VAC)	230 +/-10% (115 +/-10% under request)
Frec. (Hz)	50/60
Regulation	Microprocessor
Screen	Digital, LCD
Buttons	Membrane push buttons
Useful dimensions (mm)	365x340
Total Dimensions (mm)	455 x 410 x 445
Weight (kg)	14,7
IP	IP54

Accesories	Referencia
<i>Platform with 9 shelves for Plasma Bags - 560 x 400 mm -</i>	Under request

Ratings	
Max. Rel. humidity.	80%
Perm. Amb. T ^a	5° a 50°C
Duration of operation	100%
Fuse	1AF/250V
Protection Class	1 (protective earth)