


Slide drying

Digital



Slide Drying Units used in the preparation of microscope slides at preparation step.
 It has up to 52 slides capacity, and slides can be arranged in different orientations, hold or bending

Model	SD52-D
Reference	30000-00260
Heating	
Display	Digital
Control	Microprocesador
Nº Slides (hold/bend)	24/52
Heating Power (W)	150
T ^a range (°C)	ambT ^a +5 - 90
Resolution (°C)	0,5°C
General Data	
Useful surface (cm)	38 x 18
Power Supply	230V 50Hz*
Weight (kg)	2,3

Accessories Slide Drying	Reference
Rod	20000-0290