


RollerMix D

Digital Rotary and Oscillating Mixer with Rollers



Digital rotary mixer with rollers, controlled by microprocessor.

- Mixer with simultaneous rotative movement.
- For test tubes, bottles and closed containers.
- Speed adjustable, controlled by microprocessor .
- Four agitation rollers forming 2/3 available rails.
- Digital display LCD backlighted indicating the selected and actual speed and timer (with audible alarm).
- Membrane push-button control
- Geared synchronization with no wear.
- Automatic restart after unexpected shut down keeping the same parameters as before the incident. An incident sign appears on the display
- Quiet operation, suitable for continuous running.

Technical Specifications:

Reference	30000-00082
Agitación	
Power (W)	4
Movement	Rotación
Nº rods	4 (2 /3 rails - lines)
Speed range (rpm)	5-40
Resolution (rpm)	1
Timeer(min)	0-120
Datos Generales	
Useful length (mm)	800
Dimensions (mm)	1.000x400x110
Weight	14
IP protection	IP54